



LL 45180 S original 33x33x33
LL45190 M original 38x38x38
LL45200 L original 43x43x43

All models are available
with chest/waist strap and
crotch strap.

Lätta Lyft liftseat is soft, strong and robust. Easy to fold and take with you. The inside is glossy, which is good if you need help to move forward or back in the seat. The liftseat can also be ordered with matte inside or with both sides glossy. You do not get hot when sitting in the Lätta Lyft liftseat for long periods and you do not slide out of it. For people with spasticity, the chest strap and crotch strap can give a sense of security.



How to use Lätta Lyft liftseat

The liftseat is placed in the chair with the label facing outwards and upwards in the back. This is important in order for the lift to be correct.

Lätta Lyft's products are tested at

Technical Research Institute of Sweden, Borås, and they are CE marked.

The liftseat can hold 350 kilos, but the weight of people which can be lifted depends on whether the lift is done manually or which accessories are used. If the liftseat is connected to a patient lift, the maximum weight for the different brands will be applicable. On the liftseat there are four specially woven, easy-grip handles. It is easy to explain how lifting should be carried out. "Hold the handles." Lätta Lyft should be stretched as far as possible towards the back of the knees for the best sitting position and be straight and centred under the user. This will provide a good balance when lifting. Clothing should sit correctly when lifting and the risk for aching joints or strains because of carelessness and incorrect lifting will decrease. The front section should be lifted a little higher than the rear so the user feels secure and balanced.



There are several ways to lift

For users with weak neck muscles, it may be an advantage to lift from the front and rear. The head is then supported against the person standing behind.

In narrow passages, on the plane and on the stairs, this is the preferred technique. In other situations it is perhaps better to carry on each side. A tall and a short person lift just as well together as two people of the same size.

There can also be three or four people lifting, depending on the circumstances.

Lätta Lyft liftseats are machine-washed at 60°C. We recommend 40°C for longer life and colour retention. Stretch the straps after washing and air dry. Material: cotton/polyamide. Guarantee: 2 years from the date of the invoice.

Check that the liftseat does not have any defects. If damage does occur, the liftseat must be replaced. Manufactured in Sweden.

Svart på Vitt Reklam och Design AB / Lätta Lyft made in Sweden

Box 4067 S- 443 12 Gråbo Sweden Tel +46 707 525820 www.lattalyft.se www.facebook.com/lattalyft